




Fig. 4.4 Cycle Summary Page

The cycle summary page explains the contents of the selected cycle. Attributes such as duration, cycle brightness, and total sound groups are displayed on the right. In cycles, sounds are played concurrently and at the same volume per group. These groups play in succession, providing variation of sounds to simulate the experience of an MRI. Viewing the contents of these groups can be done using the “groups” panel, found on the left. Navigating across groups can be done using blue arrows near the top of the panel.

The user will now either cancel cycle execution using the “cancel” button, returning the user to the homepage, or continue to cycle execution by pressing the “play cycle” button in the bottom right. This will begin playing the cycle, revealing a countdown for the remaining duration of the cycle.